

FUTURE- FOCUSED INNOVATION

The Intellis™ with AdaptiveStim™ platform sets the standard for chronic pain management.



Medtronic
Further. Together

BUILT FOR A FUTURE THAT STARTS TODAY

The Intellis™ with AdaptiveStim™ platform overcomes the limitations of current spinal cord stimulation (SCS) systems with Medtronic exclusive future-focused technology, including:

Overdrive™ battery — optimized for both HD and LD energy demands


Enhanced AdaptiveStim™ — technology adjusts stimulation dose automatically

Snapshot™ reporting — based on 7 objective body positions

SureScan™ MRI — technology gives patients the most access to full body MRI

Learn more at
IntellisPainStim.com



 *Provides safe access to MRI scans anywhere on the body. Under specific conditions. Refer to product labeling for full list of conditions.

Products shown are proportionate in size to each other, but not actual size.


UNRIVALED BATTERY CHEMISTRY PERFORMANCE

Our Overdrive™ Battery Technology has a proprietary chemistry compared with SCS competitor batteries that are based on conventional Li ion chemistry. Our design allows for rapid recharge with minimal impact to the battery capacity at 9 years of use even with full discharge.

Under daily, empty to full, charge and discharge testing:

- Medtronic's Intellis Overdrive™ Battery Technology maintained 95% of the original capacity at 9 years and 3,200 cycles.^{1,2}
- The tested Competitor X spinal cord stimulators lost 30% of the capacity after only 3 years and 1,000 cycles.^{1,2}
- The Intellis devices can be recharged in ~1 hour from empty to full vs. up to 4 hours with the tested Competitor X devices.

Thus, Patients will have a more consistent recharge interval and faster recharge experience, even under the most aggressive stimulation conditions.



BATTERY PERFORMANCE

~1 hr.

recharge
from empty to full

up to
3x faster

recharge than traditional lithium
ion batteries

more than
95% battery
capacity


at 9 years, independent of therapy
parameters or recharge preferences



Patient can recharge a completely
discharged device



SNAPSHOT REPORTING DELIVERS OBJECTIVE PATIENT DATA



Integrates Medtronic AdaptiveStim™ technology to automatically adjust the therapy in response to seven unique body positions for more personalized treatment and accurate reporting

Snapshot™ transforms patient conversations from subjective to objective.

Complements self-reporting with shareable daily activity reports that validate patient usage and track their progress



All relevant therapy information, including the fluoroscopy image library and daily activity data, is stored on the Intellis™ device.



SECURE IN-DEVICE PATIENT THERAPY DATA



Provides access to data at **any location** for **real-time treatment decisions**



Integrates 128-bit encryption telemetry system to **ensure the highest levels of data security**



Allows patients the freedom to take therapy data **wherever they go**

EVOLVE WORKFLOW

Standardized guidance to simplify the trial and implant experience, and optimize patient options

- 1** **IF** lead placement spans the T9/T10 space after mapping, **THEN** consider the EvolveSM workflow deliberate dose strategy, using both HD and LD.



- 2** The suggested dose sequence begins with predefined HD settings.

HD **LD**

- 3** **DILIGENT PATIENT FOLLOW-UP** TO ASSESS FOR OPTIMAL PROGRAMMING

For more information on clinical data, visit [Medtronic.com/Evolve](https://www.medtronic.com/Evolve)

EVOLVE WORKFLOW RUNNING ON INTELLIS




Wireless external neurostimulator mimics the therapy delivered by the implanted device, improving the patient trial experience and supporting trial-to-implant conversions


Wireless programming and testing outside the sterile field simplify the intra-operative experience

Implantable neurostimulator with OverdriveTM technology improves the EvolveSM workflow experience by having faster recharge and less battery fade

FULL-BODY MRI ACCESS* WITHOUT COMPROMISE

The Intellis™ device integrates SureScan™ MRI technology for MRI access anywhere on the body — without power restrictions.

 * Under specific conditions. Requires SureScan™ MRI implantable neurostimulator and leads. Refer to product labeling for list of conditions.



SIMPLIFIES FUTURE
DIAGNOSES AND
TREATMENT WHILE
EASING PATIENT
CONCERNS










ALLOWS MRI TO
BE PERFORMED
IN NORMAL
OPERATING MODE*



GIVES PATIENTS THE
SAME UNRESTRICTED
MRI ACCESS AS
NON-IMPLANTED
PATIENTS*



Intellis™ with AdaptiveStim™ Platform

MODEL NUMBER	DESCRIPTION	PRODUCT
97715	The Intellis™ with AdaptiveStim™ Platform	
97745	Controller for patient use	
97755	Recharger for 97715 Intellis™ with AdaptiveStim™ Technology rechargeable neurostimulator	

MODEL NUMBER	DESCRIPTION	PRODUCT
THE FOLLOWING ITEMS HIGHLIGHTED IN BLUE ARE FOR THE TRIAL KIT		
97725	Wireless External Neurostimulator	
977D160	Vectris™ 1x8 subcompact trial lead kit Length: 60 cm	
977D260	Vectris™ 1x8 compact trial lead kit Length: 60 cm	
977A160	Vectris™ SureScan™ MRI 1x8 subcompact lead kit. Length: 60 cm	
977A190	Vectris™ SureScan™ MRI 1x8 subcompact lead kit. Length: 90 cm	
977A175	Vectris™ SureScan™ MRI 1x8 subcompact lead kit. Length: 75 cm	
977A260	Vectris™ SureScan™ MRI 1x8 compact. Length: 60 cm	
977A275	Vectris™ SureScan™ MRI 1x8 compact. Length: 75 cm	
977A290	Vectris™ SureScan™ MRI 1x8 compact. Length: 90 cm	

MODEL NUMBER	DESCRIPTION	PRODUCT
977C190	Specify™ SureScan™ MRI 5-6-5 surgical lead kit for spinal cord stimulation. <i>Length: 90 cm</i>	
977C165	Specify™ SureScan™ MRI 5-6-5 surgical lead kit for spinal cord stimulation. <i>Length: 65 cm</i>	
977C265	Specify™ SureScan™ MRI 2x8 surgical lead kit for spinal cord stimulation. <i>Length: 90 cm</i>	
977C290	Specify™ SureScan™ MRI 2x8 surgical lead kit for spinal cord stimulation. <i>Length: 65 cm</i>	
37500301	Wireless External Neurostimulator Boot	
CT900A	Clinician Tablet	
8880T2	Communicator	

References:

1. Remaining useful life assessment of lithium-ion batteries in implantable medical devices. Journal of Power Sources 375 (2018) 118–130.
2. Data on File.

INDICATIONS Spinal cord stimulation (SCS) is indicated as an aid in the management of chronic, intractable pain of the trunk and/or limbs-including unilateral or bilateral pain. **CONTRAINDICATIONS Diathermy** - Energy from diathermy can be transferred through the implanted system and cause tissue damage resulting in severe injury or death. **WARNINGS** Sources of electromagnetic interference (e.g., defibrillation, electrocautery, MRI, RF ablation, and therapeutic ultrasound) can interact with the system, resulting in unexpected changes in stimulation, serious patient injury or death. An implanted cardiac device (e.g., pacemaker, defibrillator) may damage a neurostimulator, and electrical pulses from the neurostimulator may cause inappropriate response of the cardiac device. **PRECAUTIONS** Safety and effectiveness has not been established for pediatric use, pregnancy, unborn fetus, or delivery. Avoid activities that put stress on the implanted neurostimulation system components. Recharging a rechargeable neurostimulator may result in skin irritation or redness near the implant site. **ADVERSE EVENTS** May include: undesirable change in stimulation (uncomfortable, jolting or shocking); hematoma, epidural hemorrhage, paralysis, seroma, infection, erosion, device malfunction or migration, pain at implant site, loss of pain relief, and other surgical risks.

Refer to www.medtronic.com for product manuals for complete indications, contraindications, warnings, precautions and potential adverse events. Rx only. Rev 0119

Medtronic

710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Tel: (763) 514-4000

medtronic.com

© 2019 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.™ *Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. UC201711476b EN